EN Mediblink Ultrasonic Nebulizer M480

INSTRUCTIONS FOR USE

Thank you very much for purchasing Mediblink Ultrasonic nebulizer M480. Be sure to read this Instruction Manual before using the unit for you to use it safely and correctly.

1

- Do not pull the power cord to unplug the adapter from the power outlet.
- Do not use extensions, unwind the power cord and keep it away from sources of heat.
- Do not tamper with the power cord for any reason. If damaged, contact your retailer.
- Always use the "U" switch to turn the device off and then unplug the adapter.
- Unplug the power cord and the adapter immediately after use and before cleaning.
- Do not immerse the device and do not pour water or other liquids onto it.
- If the device falls in the water, do not try to pick it up. Immediately unplug it. Do not reuse the device after removing it from the water.
- Only use original accessories of the manufacturer.
- Before each use, make sure that the parts have been assembled correctly.
- Do not use the device if the connecting tube is bent.
- Do not use the device close to high-frequency electromagnetic emitters.
- Keep the aeration opening free from dust and impurities. Avoid leaning the device on soft surfaces (i.e., bed or pillows).
- The device is built with materials that could cause allergic reactions (connecting tube made of EVA, PVC masks, PP mouthpiece and PET medical ampoules).
- Do not dry the device or its accessories in the microwave.

! Important notes for best operation

- Do not use the equipment, if the nebulizer water tank is without water or if the water level is below the indicated level.
- The disposable medicine cup is single use; replace it after each use to guarantee best nebulization and to avoid damage to the equipment.
- Do not use oily medicines with a viscosity over 3 Pcs (centipoise – the unit to measure viscosity), for example cortisone, as it could damage the transducer.
- Do not use essential oils (such as eucalyptus essence), active ingredients with solid elements (suspended) and viscous solutions as they could damage the transducer.
- After every use, the medicine cup and the transducer should be cleaned with a cloth dipped in pure alcohol.
- Carefully clean the upper cover of the device.
- Wait at least 15 minutes from one use to the next.


5

Table of Contents

1. Getting to Know Your Nebulizer	3
1.1 Intended Use & Indications for Use	3
1.2 Specification	3
1.3 Contraindications and Adverse reactions	4
1.4 Patient Populations	4
2. Safety Warnings	4
3. Components / Accessories	6
4. Operation Procedure	7
5. Maintenance and Disinfection	10
6. Troubleshooting Problems	11
7. Guidance and Manufacturer's Declaration	12
8. List of Symbols	14
9. Warranty	15

2

3. Components/Accessories

Name of Accessories:

- a) Child/Adult Mask
- b) Mouthpiece
- c) Adaptor
- d) Masks connector
- e) Medicine cups
- f) Connecting tube

Name of each component:

- 1) Upper cover
- 2) Medicine cup
- 3) Water tank
- 4) "U" switch
- 5) Main unit

6

1. Getting to Know Your Nebulizer

1.1 Intended Use & Indications for Use

The Mediblink Ultrasonic nebulizer M480 (model: VP-U2) is intended to aerosolize physician-prescribed solutions for inhalation (saline solution, drugs ...) that are approved for nebulization. It is not intended for life support nor does it provide any patient monitoring capabilities. Before use it is important to read the instructions and warnings contained in this booklet carefully and keep it in a safe place for future reference.

1.2 Specification

Product Name	Mediblink Ultrasonic Nebulizer M480
Model	VP-U2
Reference number	M480
Power source	AC 100~240 V, 50 Hz/60 Hz
Working Environment	Temperature: 5 to 40 °C Humidity: 15% to 90% RH, non-condensing Atmospheric Pressure: 700~1060 hPa
Storage Environment	Temperature: -25 to 70 °C Humidity: 0% to 90% RH, non-condensing Atmospheric Pressure: 700~1060 hPa
Power consumption	Adaptor: I/P: 100-240 V AC, 50/60 Hz, 0.35 W MAX.
Capacity of Medicine Cup	≤35 ml
Nebulization Rate	approx. 0.6 ml/min – 0.9 ml/min
Particle Size	MMAD approx. 3.24 μm
Percentage of particulate distribution (PDRR)	approx. 75%
Dimensions	89.9 (L) x 120 (W) x 127.4 (H) mm
Weight	Approximately 290 g (not including accessories)
Pollution	Degree 2
Overvoltage Category	CATEGORY II
Altitude	≤3000 m
Protection against harmful ingress of water or particulate matter	IP22
Ultrasonic Frequency	1.7 Mhz

3

4. Operation Procedure


1) Before each use, carefully clean and disinfect the water tank, the upper cover, the mask connector, the connecting tube, the mouthpiece and the masks according to paragraph "Maintenance and disinfection".

2) Place the device on a flat surface.

3) Remove the upper cover. Fill with tap water up to the red level indicator. This water is used as liquid conducting to the drug the ultrasonic waves and is never nebulized.


! Do not put the water tank directly under the tap to fill it with water. It is strictly forbidden to use salt solution and water with any chemical element.

Do not pour drugs directly in the water tank.

4) Fill the medicine cup with the drug, pour the drug paying attention not to exceed the maximum capacity (35 ml), then place it in the water tank. Follow the instructions for use and the dosing of the medicinal product you are using.

5) Close the water tank with the upper cover. Watch the marking on the upper cover and main unit to correctly close it.

6) Now, connect the connecting tube with the upper cover on the one end and on the other end connect it to either mouthpiece or to mask connector and mask for adults/children.

1.3 Contraindications and Adverse Reactions

Patients allergic to aerosolized drug. The device should be used with caution or be used under the guidance of doctors by those patients, who may show difficulty in breathing, apnea or continuous asthma during the process of atomization.

1.4 Patient Populations

This device is intended for use by children or adults.

2. Safety Warnings

Prior to use, ensure the product is intact and that there is no visible damage. If in doubt, contact the retailer.

- Keep the plastic bag away from children: danger of suffocation.
- This device must be used only for its intended purpose and according to the user instructions. All other types of use shall be considered inappropriate and therefore hazardous. The manufacturer shall not be deemed liable for any damage caused by inappropriate or incorrect use.
- People aged 14 younger or older, people with physical, sensory or mental disabilities or unskilled people, may use and maintain the product only under adult supervision. Children must not play with the device.
- Ensure that the operator, i.e., the patient, is able to use the appliance correctly.
- Keep all parts of the packaging as well as the connecting tube and adaptor away from children (risk of suffocation and strangling).**
- Before connecting the device, make sure the mains voltage corresponds to the voltage on plate affixed to the product.
- Do not install the appliance in areas that are difficult to access, for example too close to the wall as it would be too difficult to disconnect the device from the mains.
- Never leave the device operating unsupervised. Once finished using, turn it off and unplug it from the mains.
- Handle the product with care, protect it against accidental impact, extreme temperature fluctuations, humidity, dust, direct sunlight and sources of heat.
- In the event of a failure and/or malfunction, turn off the device without tampering with it. Always contact the retailer for repairs.
- Do not service the device during use.
- Be sure to have dry hands before connecting or disconnecting the device and when pressing the "U" switch.

4

In case of inhalation through mouthpiece:

- a) In case of troubles at low respiratory system, for example cough or bronchial infections, use the mouthpiece.
- b) Sit upright
- c) Connect the connecting tube and the upper cover, in the other end insert the mouthpiece.
- d) Tighten the mouthpiece well with the lips.
- e) Inhale slowly and deeply through the mouth and exhale through the nose. To increase the efficacy of the therapy, hold your breath for a while after having inhaled. During the treatment, remain quiet and relaxed. Breathe regularly and not too quickly. If you wish to interrupt the inhalation shortly, remove the mouthpiece from the mouth and then reinsert it.

In case of inhalation with mask:

- a) The mask for inhalation is suitable for the treatment of the nasal and pharyngeal area.
- b) Sit upright.
- c) Connect the connecting tube and the upper cover, in the other end insert the mask connector and the mask.
- d) Keep the mask lightly pressed against the nose; the mask has to cover mouth and nose without squashing them.
- e) Inhale and exhale slowly and deeply through the nose. To increase the efficacy of the therapy, hold your breath for a while after having inhaled. During the treatment, remain quiet and relaxed. Breathe regularly and not too quickly. If you wish to interrupt the inhalation shortly, remove the mask for a while and then place it again on nose and mouth.

Mouthpiece


Mask

8

7

7) Connect the adapter to the device and then to the power supply plug paying attention that the voltage is right.

8) Switch on the device pressing the "O" switch. Make sure that the warning light switches on (green color). During inhalation, the device must compulsorily remain in vertical position: this assures its correct operation, as well as the maximum efficacy of therapy. Immediately interrupt the treatment if you feel unease during the aerosol-therapy.

9) The device has an automatic stop function, which deactivates the device when the medicine is almost exhausted or after 13 minutes of inhalation. The device can be switched off at any time by pressing the "O" switch.

IMPORTANT: As with most equipment for aerosol therapy, at the end of the inhalation treatment, a certain amount of the pharmaceutical will remain in the medicine cup: this is perfectly normal. This amount of pharmaceutical, also known as residual volume, cannot be nebulized.

10) Let the device cool down for a few minutes, turn the device upside down or use an absorbing cloth to remove the residual liquid and perform the cleaning procedures according to paragraph "Cleaning and maintenance".

11) Wait at least 15 minutes from one use to the next.

5. Maintenance and Disinfection

⚠ Attention! Hot surface! Before cleaning the appliance, allow it to cool down as the transducer heats up during use.

- After removing the plug, remove the upper cover and empty any residues from the tank.
- After each use, clean the water tank and the transducer with a cloth soaked in 70% ethyl alcohol. Clean the upper cover of the appliance with a cloth soaked in water or vinegar so that no residue remains. When cleaning the upper cover, handle the flow conveyor carefully.
- Clean the outside of the appliance using a soft cloth dampened with water and mild soap (do not clean the main unit by placing it under running water and prevent any liquid from entering the appliance).
- In case of daily use, it is advisable to disinfect the accessories every day to limit the risk of infections. In this case, use a commercially available alcohol-based medical disinfectant respecting the doses and limitations provided by the disinfectant manufacturer. Immerse all the accessories in the disinfectant solution for the specified period. Completely remove the residues of the disinfectant by thoroughly rinsing the accessories under running water and allow them to air dry.
- Do not use autoclave sterilizers, EO (ethylene oxide) gas sterilizers or low temperature plasma sterilizers.
- To disinfect the appliance, it is also possible to use 70% alcohol (or an alcohol-based disinfectant spray). Check that no disinfectant residues remain in the appliance.
- The appliance, the water tank and all the accessories must be perfectly dry before using them.
- Store the appliance in a cool, dry place.
- It is advisable to check the performance of the appliance every 2 years or after a repair. Contact the Mediblink assistance service.

6. Troubleshooting Problems

Trouble	Possible Cause	How to Correct
No nebulization or low nebulization rate when the power is on.	No medication in the Medicine Cup. Too much or too little medication in the Medicine Cup.	Add the correct amount of prescribed medication to the Medicine Cup.
	The Medicine Cup is not correctly assembled.	Make sure the Medicine Cup is correctly assembled and the inhalation accessory is correctly attached.
	The device is tilted at an incorrect angle.	Do not tilt the device so the angle of the kit is greater than 45 degrees.
The device is very hot.	The connecting tube is incorrectly attached.	Make sure the connecting tube is correctly attached to the Upper cover.
	The connecting tube is folded or damaged.	Make sure the connecting tube is not folded, kinked or bent. Inspect the connecting tube for any damage. Replace the connecting tube if damaged.
The device is covered.	The device is covered.	Do not cover the device with any type of cover during use.
	Operating continuously over 13 minutes.	Limit use to 13 minutes at a time and allow a 15 minutes interval before using the device again.
No power to the device when the Power Switch is on.	The DC adapter is not plugged into the main unit.	Turn the Power Switch off. Plug the DC adapter into the DC jack. Turn the device on.

Note: If any operator requests more information such as circuit diagrams, parts list and product descriptions for repairs carried out by qualified technical personnel, please contact us.
If the unit does not nebulize normally after taking the above-mentioned procedure, contact the store where you purchased the unit or the nearest Mediblink distributor.

7. Guidance and Manufacturer's Declaration

- DEVICE needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.
- Portable and mobile RF communications equipment can affect DEVICE.
- Warning that the use of accessories, transducers and cables other than those specified with the exception of transducers and cables sold by the manufacturer of the DEVICE as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the DEVICE. Warning that the DEVICE should not be used adjacent to or stacked with other equipment. During Voltage Dips test, the product will change its operation state from on to off, after test, it needs operator to press on/off button to recover. This phenomenon does not affect BASIC SAFETY or ESSENTIAL PERFORMANCE, according to the standard, it's deemed to comply with the requirement.

A1 Electromagnetic Emissions – For all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration-electromagnetic emission		
The DEVICE is intended for use in the electromagnetic environment specified below. The customer or the user of the DEVICE should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The DEVICE uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	The DEVICE is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEC61000-3-3	Complies	

A2 Electromagnetic Immunity – For Home Healthcare Environment EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity			
The DEVICE is intended for use in the electromagnetic environment specified below. The customer or the user of the DEVICE should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to ground	±1 kV line to line ±2 kV line to ground	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U _n , 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0% U _n , 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the DEVICE requires continued operation during power mains interruptions, it is recommended that the DEVICE be powered from an uninterruptible power supply.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m	

NOTE: UT is the a.c. mains voltage prior to application of the test level.

8. List of Symbols

Symbol	Explanation
	Electrical devices are recyclable material and should not be disposed of with household waste after use! Help us protect the environment and save resources and take this device to the appropriate collection points. Please contact the organization which is responsible for waste disposal in your area if any questions.
	Applied part of type BF
	Class II appliance
	Before use, read the instructions. Electronic instructions for use: http://www.mediblink.com/fm480a.pdf
IP22	Degree of protection of electrical equipment covers, where the first figure indicates the degree of protection against penetration by solid foreign bodies (from 0 to 6) and the second figure the degree of protection against penetration by liquids (from 0 to 8).
	Manufacturer
EC REP	Authorized Representative in EU
	Warning!
LOT	Batch number*
REF	Product reference number
#	Model number
MD	Medical device
UDI	Unique Device Identifier
	Number of products included in packaging
CE	CE Marking of Conformity

*The production date can be read from LOT number [YYYYMM]; first four digits represent the Year and last two digits Month of production. Example: LOT 202503 = March 2025

9. Warranty

Product: Mediblink Ultrasonic Nebulizer M480

Manufactured for (Importer for EU and Distributor): Mediblink d.o.o., Gubčeva cesta 19, 8210 Trebnje, Slovenia; info@mediblink.com; www.mediblink.com

Seller's name, address, signature and stamp*:

Date of extradition/sales*:

*If the invoice is accompanied by this warranty, and if all above information can be seen from the invoice, it is not necessary to fill in this field.

WARRANTY TERMS

Dear customers!

The warranty period is **3 years** and starts from the day of product purchase. The warranty period on spare parts: masks, mouthpiece and tube is **1 year**. In case of product claim, you have to show the invoice. We kindly ask you to save the invoice!

Unfortunately, wrong handling with the device is a reason for 95% of customer complains. You can easily avoid any problem, by getting useful information provided by our special service department. To reach our service department, you can call or send e-mail to Mediblink local Distributor.

Before sending the product back to retailer, we kindly ask you to call our service department, to get help about how to use the device to save you with unneeded trips.

The Manufacturer guarantees free elimination of all imperfections due to defects in material or manufacturing procedure by repairing or replacing the product. In case that the product cannot be repaired or replaced, the customer will get the money refund. The guarantee is not valid in case of the force majeure, accidents or unexpected events (such as lightning, water, fire etc.), incorrect use or incorrect transport, non-compliance with safety and maintaining regulations or in case of unprofessional product intervention.

Traces of every day product usage (scratches, abrasions) and not subject to claim. The warranty does not eliminate the customer rights, which originate from seller responsibility for product flaws. By accepting the claimed product by the service department, the service department does not take responsibility for loss of saved data or settings on the product. All product repairs, which are performed out of product warranty period, have to be paid by customer by prior notice.

The Manufacturer guarantees the product quality and flawless product operation in the warranty period, which starts with the day of product purchase. If the product cannot be repaired in 45 days, the product will be replaced with a new one. In case that the product cannot be replaced, the money will be refunded to the customer.

In case of product claim, call or send e-mail to Mediblink local Distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/ or patient is established.

Replacement parts can be ordered at www.mediblink.com.

Manufacturer:
Vapo Healthcare Co., Ltd.
Southern unit of third floor,
building B, No. 99 Yudai West Rd,
High tech district, Kunshan Suzhou,
215301 Jiangsu, P.R. China

Manufactured for (Importer for EU & Distributor):
Mediblink d.o.o., Gubčeva cesta 19,
8210 Trebnje, Slovenia;
info@mediblink.com;
www.mediblink.com

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